
Characteristics / Caratteristiche

Working surface / Piano lavoro

- The tissue embedding cassettes and base molds trays are removable and can be set at different temperatures
- The heated forceps holder can be set at the preferred temperature.
- The ergonomic large working surface allows processing many samples at the same time.
- The setting knobs allow to adjust the precise paraffin flow.

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

- The large cold spot Peltier is designed to orient even small and big samples (Mega Cassettes).
- The temperature range is from room temp. To -5 °C.
- The large work surface allows draining of the excess paraffin into two large collecting trays and allows processing many samples at the same time.
- The temperature can be set from ambient to +99 °C and can be adjusted through the control panel.

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All functions of the TEC 2800 are settable through the control panel

- Temperature of working surface
- Temperature of molds cassettes warming trays
- Temperature of paraffin tank
- Work time and day can be programmed automatically

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Cryo Console

The TEC 2800 Cryo Console is composed of a compressor cooled work surface held to chill the paraffin blocks rapidly. The Teflon coated surface allows a minimum ice formation and easy cleaning. Temperature is adjustable from ambient to -13 °C. TEC 2800 Cryo Console can be used as a stand alone unit near the microtome station.

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Thermal Console

The TEC 2800 Thermal Console holds processed tissue specimens, in any type of cassettes, awaiting paraffin embedment. Stores and warms samples in a large paraffin tank "2,5 Lt." Temperature can be set from ambient to 99 °C.

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Histo-Trim

The TEC 2800 Histo Trim is a practical device designed to remove excess paraffin from the embedding cassettes after inclusion of specimens. Scraping the embedded block on the "grooved plane" heated surface, excess paraffin will be removed leaving the cassette clean ready to be set properly on the clamp of the microtome. Temperature can be set from ambient to 99 °C.

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Foot Switch / Pedale interruttore

TEC 2800

Nominal supply voltage 220 V/50 Hz
Maximum power draw 1000 VA
Working temperature from ambient to 70 °C
increments 1 °C

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Cryo Console

Nominal supply voltage 220 V/50 Hz
Maximum power draw 110 VA
Working temperature from ambient to -13 °C
increments 1 °C

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Thermal Console

Nominal supply voltage 220 V/ 50 Hz
Maximum power draw 400 VA
Working temperature from ambient to 70 °C
increments 1 °C

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