



MODEL GL100 ELEVATING PATHOLOGY WORK STATION

MODEL GL100

Model GL100 Elevating Pathology Gross Station comes standard with all required features for any grossing procedure and has a list of optional accessories to suit anyone's personal needs. The unit features all stainless steel construction of the entire work area and all exterior closure panels.

Standard Design Features

- All stainless steel construction
- Height adjustment from 32 to 44 via electro-hydraulic mechanism
- Exhaust duct for connection to building ventilation system with adjustable bellows or provided as self-contained with internal exhaust blowers and filter system
- Units features both downdraft and backdraft ventilation
- Foot operated faucet control for hot and cold water
- Spray hose assembly with hand control
- 1/2 hp commercial disposal
- Dissecting area rinse to provide a constant flow of water to work area
- Top mounted fluorescent light fixture
- Stainless C-fold paper towel holder
- Magnetic instrument holder
- Polyethylene dissecting board
- GFCI waterproof electrical receptacle





MODEL GL100 LAYOUT DRAWING, OPTIONAL FEATURES

MODEL GL100

Optional Accessories:

- Heavy-Duty commercial disposal (3/4 hp, 1 hp, 2 hp)
- Self-Contained ventilation system with replaceable filters
- Gas / Air / Vacuum Valves
- Digital platform scale, 2 kg capacity, economy or stainless grades available
- Dictation Equipment with foot control or voice operation
- Camera mount facilities
- Eyewash assembly
- Cassette holders
- Computer keyboard and monitor stand
- Plexiglass splash shield to enclose work area
- Forms holder
- Formalin dispensing system for new and waste formalin
- Removable stainless steel ruler
- Magnifier light on flex arm
- Halogen light on flex arm

