



### MODEL GL106 ELEVATING COUNTER GROSS STATION



### MODEL GL106

Model GL 106 Elevating Counter Gross Station comes standard with all required features for any grossing procedure and has a list of optional accessories to suit anyone's personal needs. The unit features all stainless steel construction of the entire work area and all exterior closure panels.

#### Standard Design Features:

- All stainless steel construction
- Height adjustment from 32" to 44" via electro-mechanical actuator
- Exhaust duct for connection to building ventilation system
- Large integral sink with mixing faucet
- Foot operated control of faucet for hot and cold water
- Spray hose assembly with hand control
- 1/2 hp commercial disposal
- Small specimen rinse basket with stainless steel holder
- Incandescent light with 3x magnifier mounted on a flexible arm
- Formalin dispensing and collection system
- Magnetic instrument holder
- Polyethylene dissecting board
- GFCI waterproof electrical receptacles



**MODEL GL106 LAYOUT DRAWING, OPTIONAL FEATURES**

**MODEL GL106**

**Optional Accessories:**

- Heavy-Duty commercial disposal (3/4 hp, 1 hp)
- Gas / Air Valves / Vacuum / Oxygen
- Eyewash Assembly

