



MODEL GL110 COUNTERTOP PATHOLOGY WORK STATION



MODEL GL110

An economical bench pathology workstation with all the necessary fixtures and accessories to perform routine grossing procedures.

Standard Design Features:

- Available in two sizes: 48 inches and 60 inches long
- All stainless steel construction
- Integral sink with mixing faucet
- Top mounted fluorescent light fixture
- Two-incandescent lights with 3x magnifier mounted on flexible arms
- Formalin container with spigot
- C-fold paper towel holder, stainless steel
- Magnetic instrument holder
- Polyethylene dissection board

**MODEL GL110 LAYOUT DRAWING, OPTIONAL FEATURES & FILTERS****MODEL GL110****Optional Accessories:**

- Self-Contained Ventilation system with replaceable filters
- Digital Platform Scale, 2 kg capacity, economy or stainless grades available
- Dictation Equipment with foot control or voice operation
- Camera mount facilities
- Eyewash Assembly
- Model GL120 or GL125 Station Bases

Replacement Filters:

- Model GL 112 Activated Charcoal Filter for odorous fumes
- Model GL 113 Potassium Permanganate for formaldehyde fumes

