

**MODEL GL117 PATHOLOGY WORK STATION WITH REAR EXHAUST****MODEL GL117**
Standard Design Features

- Sink top, basin, drainboards and backsplash made from welded, 14 gauge stainless steel with a No. 4 finish
- Units can be fabricated in varying lengths, depths and heights
- Specimen shelf, welded heavy gauge stainless steel
- Removable 16" x 31" stainless steel dissecting grid plate, located over the drain trough
- Flushing system provided under the dissecting grid plate
- Heavy-Duty Commercial Waste Disposal (1/2 hp) with vacuum breaker and solenoid
- Spray hose assembly
- Vacuum Breaker protection to prevent backflow of contaminated water
- Hot and Cold Water Mixing Faucet with Gooseneck Spout and Wrist Blade Handles
- Fluorescent, vapor-proof light with switch
- Dissecting area is sloped toward sink basin for rapid and positive drainage
- Heavy-gauge polished stainless steel leg assembly with welded cross rails to handle heavy-duty usage
~wheel chair accessible

