



MODEL GL119 SIDE LATERAL EXHAUST PATHOLOGY WORK STATION



MODEL GL119

The Model GL 119 Pathology Work Station provides side lateral downward exhaust to the dissecting area. The unit can be operated from either side.

Standard Design Features

- All stainless steel construction
- Exhaust duct for connection to building system
- Downward exhaust on both sides of dissecting area
- Integral sink with mixing faucet
- Hot and cold hand spray assembly
- 1/2 hp commercial disposal with waterproof switch
- Dissecting area rinse assembly
- Raised perforated dissecting trays
- (4) Instrument drawers
- Pull out writing platform
- Polyethylene dissecting board
- GFCI waterproof duplex receptacle





MODEL GL119 LAYOUT DRAWING, OPTIONAL FEATURES

MODEL GL119
Optional Accessories

- Heavy-Duty commercial disposal (3/4 hp, 1 hp, 2 hp)
- Gas / Air Valves / Vacuum / Oxygen
- Eyewash Assembly
- Deck Mounted Magnifier Light

